

Disposable Hot Biopsy Forceps

Indications

In endoscopic surgery, high-frequency current is used to excise and collect tissue samples in the digestive tract, and electrocoagulation to stop bleeding.



Features

- ※ Rich design of forceps cups meets various needs of doctors
- ※ The head of the biopsy forceps can be opened and closed easily and flexibly, and the incision is sharp. Used with high frequency generator, it has the function of coagulation and electric cutting
- ※ The outer surface is coated with a special material super-smooth coating, which can be inserted smoothly into the forceps channel, effectively reducing the wear of the instruments on the endoscope forceps channel.
- ※ Small wound surface, no tissue charring, fast healing
- ※ Unique four-bar plier head design, the jaws are firmly engaged and not easy to slip off
- ※ Temperature controllable, suitable for a variety of cuts, fast coagulation, short operation time
- ※ Sterile packaging, easy to use

Specifications

Code	Jaw O.D.(mm)	Working Length (mm)	Min. access channel (mm)
HF18-24/60	2.4±0.2	600±100	2.8
HF18-24/120	2.4±0.2	1200±100	2.8
HF18-24/160	2.4±0.2	1600±100	2.8
HF18-24/180	2.4±0.2	1800±100	2.8
HF18-24/200	2.4±0.2	2000±100	2.8
HF18-24/230	2.4±0.2	2300±100	2.8
HF18-18/60	1.8±0.2	600±100	2.0
HF18-18/120	1.8±0.2	1200±100	2.0
HF18-18/160	1.8±0.2	1600±100	2.0
HF18-18/180	1.8±0.2	1800±100	2.0
HF18-18/200	1.8±0.2	2000±100	2.0